





Using Fentanyl Test Strips

1. Add water to the empty baggie or cooker (after drawing up your shot). For pills, crumble a small amount and add water
2. Mix the water well
3. Dip test strip in the water for 15 seconds
4. Put test strip on top of cooker or another clean surface
5. Wait about 5 minutes
6. Read test strip
7. If your drugs have fentanyl in them, take precautions. If you are unsure about your results, test again



Attention: If you are testing meth

Add ½ cup water to the meth residue (no shards), then run the test. This will give you the most accurate results.



This document was developed [in part] under grant number TI083243 from the Substance Abuse and Mental Health Services Administration (SAMHSA), U.S. Department of Health and Human Services (HHS). The views, policies, and opinions expressed are those of the authors and do not necessarily reflect those of SAMHSA or HHS.



**Tribal Opioid
RESPONSE**
Healing our Nations Together



NPAIHB
Indian Leadership for Indian Health

www.npaihb.org/opioid